

## Find Host ID (MAC/ physical address) for IQTL license activation

### Find Host ID for Windows computers

STEP	ACTION
1	In the Windows search bar, type <code>cmd</code> and press <b>Enter</b> on your keyboard
2	Type <code>ipconfig/all</code> and press <b>Enter</b>
3	Select a network adapter for the current system: <ul style="list-style-type: none"><li>• Via an Ethernet cable, choose <b>Ethernet Adapter Ethernet</b></li><li>• Via a wireless network adapter, choose <b>Wireless LAN adapter Wi-Fi</b></li></ul>
4	Find the <b>Host ID</b> as a string of numbers and letters under <b>Physical Address</b>

### Find Host ID for MAC computers

STEP	ACTION
1	Open the Apple menu and click the Apple icon
2	Scroll to <b>System Preferences</b>
3	Select Network → Wi-Fi → Advanced
4	Copy the <b>MAC</b> address shown under the <b>Wi-Fi Address</b> or <b>Airport Address</b>